

DIGITAL INCUBATOR



Construction : Incubator is made up of double walled construction. Front door with gasket and lock. The motorized blower is fitted on the backside for uniform temperature (only forced air circulation model). Insulation is provided by means of 70mm. Glass Wool insulation. The unit is provided with adjustable stainless steel rod trays.

Also Available in

- ◆ **Customized Sizes**
- ◆ **GMP model**

MODEL SPECIFICATION :

MODEL	SIZE (mm) W x D x H	LITRE	RATING (Watts)
PINS -350	350 x 350 x 350	42	1000
PINS -450	450 x 450 x 450	90	1500
PINS -600	600 x 600 x 600	216	2000

TECHNICAL SPECIFICATION:

Temp. Controller	PID Controller
Temp. Sensor	PT100 - RTD
Temp. Range	+5°C above ambient upto 70°C
Temperature accuracy	±0.5°C
Inner Chamber	Made of Stainless steel Grade - 304
Outer Chamber	Made of Mild steel with powder coated
Insulation	70mm Glass wool Insulation
Supply	230V AC, Single Phase
Blower motor	1/70HP AC, 1340 RPM