



FUME HOOD

Construction : Made of mild steel body with epoxy powder coated FRP COATING INNER CHAMBER With working table covered of stainless steel 316 or 304 quality. Acid resistant. WORKING PLATFORM ON 18MM THICKNESS BLACK GRANITE. The work surface is fitted with small wash basin of size 250 x 250 x 125mm depth of stainless steel.(optional). A gas cock is also fitted for gas supply. Chamber is provided with fluorescent lighting Socket with switches for light and exhaust system Front opening of the fume hood will be H 585 x W 910mm.

Also Available in

- ◆ Customized Sizes
- ◆ GMP model



MODEL SPECIFICATION :

MODEL	SIZE (MM) W X D X H	RATING(WATTS)
PFH-2	600 x 600 x 1500 mm	300
PFH-3	900 x 600 x 1500 mm	600
PFH-4	1200 x 600 x 1500 mm	750
PFH-5	1500 x 600 x 1500 mm	1200