



WATER BATH

Construction : Made Of Mild Steel Body And Inner Made Of Stainless Steel Sheet 304 Quality. Concentric Rings Are Also Made Of Stainless Steel & Polymide Rings. Heating Elements Have An Automatic Cut Off Device When The Bath Gets Dry. With Arrangement Of 6 Holes Or 12 Holes. Electrically Operated On 220/230v AC.

Also Available in

- ◆ Customized Sizes
- ◆ GMP model



MODEL SPECIFICATION :

MODEL	SIZE	HOLES DIA	LITRE	RATING
PWB-6	250 x 350 x 100	75 DIA	8	750 W
PWB-12 1000 W	300 x 400 x 100	75 DIA	12	1000 w

TECHNICAL SPECIFICATION:

Temp.Controller	PID Controller
Temperature range	5°C to 95°C
Temperature accuracy	±0.5°C
Temp.Sensor	PT100 - RTD
Inner Chamber	Made of Stainless steel Grade - 304
Outer Chamber	Made of Mild steel with powder coated