



P-LAB

LAB ESSENTIALS

For Every Laboratory



PRECISION LAB FURNITURE & INSTRUMENT



HOT AIR OVEN

Construction : Hot air oven is made up of tripled walled construction. The cabinets are made using stainless steel material.

The 70mm thick glass wool insulation is provided. The heat is generated inside the oven using the heater fitted on inside the chamber.

Also Available in

- ◆ **Customized Sizes**
- ◆ **GMP model**



MODEL SPECIFICATION :

MODEL	SIZE (mm) W x D x H	LITRE	RATING (Watts)
PHOS -350	350 x 350 x 350	42	1000
PHOS -450	450 x 450 x 450	90	1500
PHOS -600	600 x 600 x 600	216	2000

TECHNICAL SPECIFICATION:

Temp. Controller	PID Controller
Temp. Sensor	J Type Thermocouple
Temp. Range	+ 5°C above ambient upto 225°C
Temperature accuracy	±1°C
Inner Chamber	Made of Stainless steel Grade - 304
Outer Chamber	Made of Mild steel with powder coated
Insulation	70mm Glass wool Insulation
Supply	230V AC, Single Phase
Blower motor	1/70HP AC, 1340 RPM